

Anti-Exoc2 Rabbit Monoclonal Antibody For Wb, Elisa - Q96Kp1

Item Number: CM0007095

Introduction

Recombinant rabbit monoclonal antibody targeting human EXOC2 (exocyst complex component 2, Sec5). Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying exocytosis and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	EXOC2 Rabbit mAb
Remarks/Alias	SEC5; Sec5p; SEC5L1; NEDFACH; EXOC2
Species	Human
Gene ID (Human)	55770
Gene ID	55770
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KP1
Protein Weight	104 kDa
Shipping	Ice bag