

Anti-Exoc2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q96Kp1

Item Number: CM0023926

Introduction

Rabbit polyclonal antibody targeting EXOC2/Sec5, validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for exocyst complex and vesicle trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	EXOC2 Rabbit pAb
Remarks/Alias	SEC5; Sec5p; SEC5L1; NEDFACH; EXOC2
Species	Human
GeneID (Human)	55770
GeneID	55770
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KP1
Protein Weight	104kDa
Shipping	Ice bag