

Anti-Exoc3 Monoclonal Antibody For Wb, Elisa - O60645

Item Number: CM0003251

Introduction

EXOC3 Rabbit mAb (O60645) is a validated monoclonal antibody for WB, ELISA with human, mouse, rat reactivity. Investigate exocyst complex-mediated vesicle docking and secretion in cellular transport.

[Learn More](#)

Parameter	Value
Product Name	EXOC3 Rabbit mAb
Remarks/Alias	SEC6; Sec6p; SEC6L1; EXOC3
Species	Human
GeneID	11336 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 650-745 of human EXOC3 (O60645)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60645
Protein Weight	86kDa
Shipping	Ice bag