

Anti-Exoc3 Polyclonal Antibody For Wb, Elisa - O60645

Item Number: CM0027682

Introduction

Rabbit polyclonal antibody targeting human EXOC3 (Sec6), a component of the exocyst complex. Suitable for Western blot and ELISA applications in exocytosis and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	EXOC3 Rabbit pAb
Remarks/Alias	SEC6; Sec6p; SEC6L1; EXOC3
Species	Human
Gene ID (Human)	11336
Gene ID	11336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 396-745 of human EXOC3 (NP_009208.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60645
Protein Weight	86kDa
Shipping	Ice bag