

Anti-Exoc4 Polyclonal Antibody For Wb And Elisa - Q96A65

Item Number: CM0013048



Introduction

Anti-EXOC4 rabbit polyclonal antibody validated for Western blot and ELISA. Targets exocyst complex component Sec8 in human, mouse, and rat. Ideal for exocytosis and protein transport research.

[Learn More](#)

Parameter	Value
Product Name	EXOC4 Rabbit pAb
Remarks/Alias	SEC8; Sec8p; SEC8L1; EXOC4
Species	Human
Gene ID (Human)	60412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human EXOC4 (NP_068579.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96A65
Protein Weight	110kDa
Shipping	Ice bag