

Anti-Exoc5 Rabbit Polyclonal Antibody For Wb And Elisa - O00471

Item Number: CM0029421

Introduction

This rabbit polyclonal antibody targets human EXOC5, a key exocyst complex component. It is validated for WB and ELISA and cross-reacts with mouse and rat, making it suitable for studying exocytic vesicle docking and secretion.

[Learn More](#)

Parameter	Value
Product Name	EXOC5 Rabbit pAb
Remarks/Alias	SEC10, HSEC10, SEC10P, PRO1912, SEC10L1, EXOC5
Species	Human
Gene ID (Human)	10640
Gene ID	10640
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-610 of human EXOC5 (NP_006535.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00471
Protein Weight	82kDa
Shipping	Ice bag