

Anti-Exosc1 Polyclonal Antibody For Wb And Elisa - Q9Y3B2

Item Number: CM0013809



Introduction

EXOSC1 Rabbit Polyclonal Antibody (Q9Y3B2) for WB and ELISA detects human exosome component CSL4, a non-catalytic subunit of the RNA exosome complex involved in RNA processing and degradation.

[Learn More](#)

Parameter	Value
Product Name	EXOSC1 Rabbit pAb
Remarks/Alias	p13; CSL4; SKI4; Csl4p; PCH1F; Ski4p; CGI-108; EXOSC1
Species	Human
Gene ID (Human)	51013
Gene ID	51013
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human EXOSC1 (NP_057130.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y3B2
Protein Weight	21kDa
Shipping	Ice bag