


Anti-Exosc10 Polyclonal Antibody For Wb And Elisa - Q01780

Item Number: CM0023895



Introduction

Rabbit polyclonal antibody against human EXOSC10, validated for WB and ELISA. Targets the RNA exosome component involved in RNA processing, DNA repair, and immune cell function. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	EXOSC10 Rabbit pAb
Remarks/Alias	p2; p3; p4; RRP6; PMSCL; Rrp6p; PM-Scl; PMSCL2; PM/Scl-100; EXOSC10
Species	Human
GeneID	5394
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 586-885 of human EXOSC10 (NP_001001998.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01780
Protein Weight	101kDa
Shipping	Ice bag