

Anti-Exosc3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqt5

Item Number: CM0013923



Introduction

High-quality rabbit polyclonal antibody targeting EXOSC3 (RRP40/p10), an RNA exosome core component. Validated for Western blotting and ELISA, with cross-reactivity in human, mouse, and rat samples. Ideal for RNA metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	EXOSC3 Rabbit pAb
Remarks/Alias	p10; PCH1B; RRP40; Rrp40p; CGI-102; hRrp-40; bA3J10.7; EXOSC3
Species	Human
GeneID (Human)	51010
GeneID	51010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human EXOSC3 (NP_057126.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQT5
Protein Weight	30kDa
Shipping	Ice bag