

Anti-Exosc7 Polyclonal Antibody For Wb, Elisa - Q15024

Item Number: CM0027695



Introduction

Rabbit polyclonal antibody against human EXOSC7 (RRP42), for WB and ELISA applications. Validated in human, mouse, and rat samples. Target is a non-catalytic component of the RNA exosome complex involved in RNA processing and degradation.

[Learn More](#)

Parameter	Value
Product Name	EXOSC7 Rabbit pAb
Remarks/Alias	p8; EAP1; RRP42; Rrp42p; hRrp42p; EXOSC7
Species	Human
Gene ID (Human)	23016
Gene ID	23016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-291 of human EXOSC7 (NP_055819.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15024
Protein Weight	32kDa
Shipping	Ice bag