

Anti-Exosc8 Rabbit Polyclonal Antibody For Wb, Elisa - Q96B26

Item Number: CM0030470



Introduction

Rabbit anti-EXOSC8 polyclonal antibody validated for WB and ELISA. Targets exosome component RRP43, involved in RNA processing and degradation. Ideal for human sample research.

[Learn More](#)

Parameter	Value
Product Name	EXOSC8 Rabbit pAb
Remarks/Alias	p9; CIP3; EAP2; OIP2; PCH1C; RRP43; Rrp43p; bA421P11.3; EXOSC8
Species	Human
Gene ID (Human)	11340
Gene ID	11340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-276 of human EXOSC8 (NP_852480.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96B26
Protein Weight	30kDa
Shipping	Ice bag