

Anti-Exportin 2 (Xpo2) Rabbit Monoclonal Antibody For Wb And Elisa - P55060

Item Number: CM0002332

Introduction

Rabbit monoclonal antibody against Exportin 2 (XPO2/CSE1L), suitable for Western blot and ELISA applications. This antibody detects human and mouse Exportin 2 protein (110 kDa) and is an essential tool for investigating nucleocytoplasmic transport mechanisms.

[Learn More](#)

Parameter	Description
Product Name	Exportin 2 (XPO2) Rabbit mAb
Remarks/Alias	CAS; CSE1; XPO2; Exportin 2 (XPO2)
Species	Human
Gene ID (Human)	1434
Gene ID	1434
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 872-971 of human Exportin 2 (XPO2) (P55060)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55060
Protein Weight	110kDa
Shipping	Ice bag