

Anti-Exportin 5 (Xpo5) Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Hav4

Item Number: CM0028394



Introduction

Rabbit polyclonal antibody to Exportin 5 (XPO5) suitable for WB, IF/ICC, IP, ELISA. Validated in human, mouse, rat. 136kDa, immunogen aa 925-1204. Critical for miRNA and dsRNA nuclear export.

[Learn More](#)

Parameter	Value
Product Name	Exportin 5 (XPO5) Rabbit pAb
Remarks/Alias	exp5; Exportin 5 (XPO5)
Species	Human
GeneID (human)	57510
GeneID	57510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 925-1204 of human Exportin 5 (NP_065801.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAV4
Protein Weight	136kDa
Shipping	Ice bag