

Anti-Exportin 5 (Xpo5) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Hav4

Item Number: CM0011658



Introduction

Rabbit monoclonal antibody against human Exportin 5 (XPO5). Validated for WB, IHC-P, ELISA and cross-reacts with mouse and rat. Ideal for studying nuclear export and miRNA biology.

[Learn More](#)

Parameter	Specification
Product Name	Exportin 5 (XPO5) Rabbit mAb
Remarks/Alias	exp5; Exportin 5 (XPO5)
Species	Human
GeneID (Human)	57510
GeneID	57510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Exportin 5 (XPO5) (Q9HAV4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAV4
Protein Weight	136kDa
Shipping	Ice bag