

Anti-Exportin-6 (Xpo6) Rabbit Polyclonal Antibody For Wb And Elisa - Q96Qu8

Item Number: CM0015690

Introduction

XPO6 Rabbit pAb detects Exportin-6 in human and mouse samples via WB and ELISA. Ideal for nuclear transport and actin dynamics studies. Polyclonal antibody raised against recombinant human XPO6 fragment.

[Learn More](#)

Parameter	Specification
Product Name	XPO6 Rabbit pAb
Remarks/Alias	EXP6; RANBP20; XPO6
Species	Human
GenelD (Human)	23214
GenelD	23214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-530 of human XPO6 (NP_055986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96QU8
Protein Weight	129kDa
Shipping	Ice bag