

Anti-Ext1 Rabbit Polyclonal Antibody For Wb, Elisa - Q16394

Item Number: CM0021585



Introduction

Rabbit polyclonal antibody against human exostosin-1 (EXT1), validated for WB and ELISA, with cross-reactivity to rat. Ideal for studying heparan sulfate biosynthesis and hereditary multiple exostoses.

[Learn More](#)

Parameter	Value
Product Name	EXT1 Rabbit pAb
Remarks/Alias	EXT; LGS; TTV; LGCR; TRPS2; EXT1
Species	Human
Gene ID (Human)	2131
Gene ID	2131
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human EXT1 (NP_000118.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q16394
Protein Weight	86kDa
Shipping	Ice bag