

Anti-Ext1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16394

Item Number: CM0020396



Introduction

EXT1 Rabbit Polyclonal Antibody (pAb) for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets exostosin-1, a heparan sulfate co-polymerase subunit critical for glycan elongation and extracellular matrix homeostasis.

[Learn More](#)

Parameter	Value
Product Name	EXT1 Rabbit pAb
Remarks/Alias	EXT; LGS; TTV; LGCR; TRPS2; EXT1
Species	Human
GeneID (Human)	2131
GeneID	2131
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 467-746 of human EXT1 (NP_000118.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16394
Protein Weight	86kDa
Shipping	Ice bag