

Anti-Ext2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q93063

Item Number: CM0022382



Introduction

Rabbit polyclonal antibody targeting human EXT2 (exostosin-2), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for HS biosynthesis and multiple exostoses research.

[Learn More](#)

Parameter	Value
Product Name	EXT2 Rabbit pAb
Remarks/Alias	SOTV; SSMS; EXT2
Species	Human
Gene ID	2132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-340 of human EXT2 (NP_000392.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93063
Protein Weight	82kDa
Shipping	Ice bag