

Anti-Eya1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99502

Item Number: CM0030530

Introduction

Rabbit polyclonal antibody targeting EYA1, a transcription coactivator and phosphatase involved in DNA repair and organogenesis. Validated for WB, IHC-P, ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	EYA1 Rabbit pAb
Remarks/Alias	BOP; BOR; BOS1; OFC1; EYA1
Species	Human
Gene ID (Human)	2138
Gene ID	2138
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human EYA1 (NP_000494.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99502
Protein Weight	65kDa
Shipping	Ice bag