


Anti-Eya3 Rabbit Polyclonal Antibody For Wb And Elisa - Q99504

Item Number: CM0030518



Introduction

Anti-EYA3 Rabbit Polyclonal Antibody validated for WB and ELISA; reacts with human, mouse, rat. Targets protein phosphatase EYA3, a key regulator in DNA repair and transcription. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	EYA3 Rabbit pAb
Remarks/Alias	EYA3
Species	Human
GeneID (human)	2140
GeneID	2140
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EYA3 (NP_001981.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99504
Protein Weight	63kDa
Shipping	Ice bag