


Anti-Eya3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99504

Item Number: CM0017712



Introduction

High-specificity rabbit polyclonal antibody targeting human EYA3 phosphatase (UniProt Q99504). Validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Confirmed cross-reactivity with rat. Ideal for studying DNA damage repair, transcription regulation, and developmental biology pathways.

[Learn More](#)

Parameter	Value
Product Name	EYA3 Rabbit pAb
Remarks/Alias	EYA3
Species	Human
Gene ID (Human)	2140
Gene ID	2140
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EYA3 (NP_001981.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q99504
Protein Weight	63kDa
Shipping	Ice bag