

Anti-Ezh1 Polyclonal Antibody For Wb And Elisa - Q92800

Item Number: CM0028250



Introduction

EZH1 Rabbit Polyclonal Antibody validated for WB and ELISA, detecting human, mouse, rat EZH1. This antibody targets a histone methyltransferase critical for transcriptional repression and embryonic stem cell maintenance.

[Learn More](#)

Parameter	Value
Product Name	EZH1 Rabbit pAb
Remarks/Alias	KMT6B; EZH1
Species	Human
GeneID (Human)	2145
GeneID	2145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-280 of human EZH1 (NP_001982.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92800
Protein Weight	85kDa
Shipping	Ice bag