


## Anti-Ezh2/Kmt6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15910

Item Number: CM0016830



### Introduction

Polyclonal antibody targeting human EZH2/KMT6 (aa 1-250). Validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Suitable for epigenetic and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EZH2/KMT6 Rabbit pAb                                                                                                     |
| Remarks/Alias    | WVS; ENX1; KMT6; WVS2; ENX-1; EZH2b; KMT6A; EZH2/KMT6                                                                    |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 2146                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EZH2/KMT6 (NP_001190176.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                             |
| SWISS            | Q15910                                                                                                                   |
| Protein Weight   | 85kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |