


## Anti-Ezh2/Kmt6 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q15910

Item Number: CM0019659



### Introduction

High-quality rabbit polyclonal antibody targeting EZH2/KMT6 (UniProt Q15910), suitable for western blotting, immunoprecipitation, and ELISA. Detects human and mouse EZH2, a key histone methyltransferase involved in transcriptional repression and cancer.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EZH2/KMT6 Rabbit pAb                                                                                                                     |
| Remarks/Alias    | WVS; ENX1; KMT6; WVS2; ENX-1; EZH2b; KMT6A; EZH2/KMT6                                                                                    |
| Species          | Human                                                                                                                                    |
| GeneID (human)   | 2146                                                                                                                                     |
| GeneID           | 2146                                                                                                                                     |
| Immunogen        | Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 158-261 of EZH2/KMT6 (NP_001190176.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                                    |
| Application      | WB, IP, ELISA                                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                                             |
| SWISS            | Q15910                                                                                                                                   |
| Protein Weight   | 85kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |