


Anti-Ezrin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P15311

Item Number: CM0003103



Introduction

High-quality Anti-Ezrin Rabbit Monoclonal Antibody for Western blot, IHC-P, and ELISA. Validated in human samples, targets Ezrin (P15311), a cytoskeletal linker protein. Ideal for research in epithelial cell biology.

[Learn More](#)

Parameter	Value
Product Name	Ezrin Rabbit mAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; Ezrin
Species	Human
GeneID (Human)	7430
GeneID	7430
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 300-400 of human Ezrin (P15311).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P15311
Protein Weight	69kDa
Shipping	Ice bag