


Anti-Ezrin/Radixin/Moesin Polyclonal Antibody For Wb And Elisa - P15311 / P35241 / P26038

Item Number: CM0019318



Introduction

High-quality rabbit polyclonal antibody targeting Ezrin, Radixin, and Moesin (ERM family). Validated for Western blot (WB) and ELISA; cross-reacts with human, mouse, and rat. Suitable for membrane-cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Ezrin/Radixin/Moesin Rabbit pAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; DFNB24; HEL70; IMD50; Ezrin/Radixin/Moesin
Species	Human
GeneID (Human)	7430/5962/4478
GeneID	7430/5962/4478
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 487-586 of human Ezrin/Radixin/Moesin (NP_003370.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15311/P35241/P26038
Protein Weight	23-71 kDa
Shipping	Ice bag