


Anti-F13A1 Polyclonal Antibody For Wb, If-P, Elisa - P00488

Item Number: CM0012896



Introduction

Rabbit polyclonal antibody against coagulation factor XIII A chain (F13A1) for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Stabilizes fibrin clots via transglutaminase activity.

[Learn More](#)

Parameter	Value
Product Name	F13A1 Rabbit pAb
Remarks/Alias	F13A; F13A1
Species	Human
GeneID (Human)	2162
GeneID	2162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-732 of human F13A1 (NP_000120.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00488
Protein Weight	83kDa
Shipping	Ice bag