

Anti-F2RL3 Polyclonal Antibody For Wb, Elisa - Q96Ri0

Item Number: CM0029989



Introduction

Rabbit polyclonal antibody against human Proteinase-activated receptor 4 (F2RL3/PA4) for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat samples. Suitable for thrombin signaling and platelet activation research.

[Learn More](#)

Parameter	Value
Product Name	F2RL3 Rabbit pAb
Remarks/Alias	PAR4; F2RL3
Species	Human
Gene ID (Human)	9002
Gene ID	9002
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-150 of human F2RL3 (NP_003941.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RI0
Protein Weight	41kDa
Shipping	Ice bag