

Anti-F4/80 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q61549

Item Number: CM0005472

Introduction

F4/80 Rabbit monoclonal antibody for use in WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Supports macrophage biology, cell adhesion, and regulatory T-cell development studies.

[Learn More](#)

Parameter	Value
Product Name	F4/80 Rabbit mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
GenelD (Human)	2015
GenelD (Mouse)	13733
Immunogen	Recombinant protein of mouse F4/80
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag