

Anti-Fa2H Rabbit Polyclonal Antibody For Wb, Elisa - Q7L5A8

Item Number: CM0015859

Introduction

Rabbit polyclonal antibody targeting human fatty acid 2-hydroxylase (FA2H). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Supports research in sphingolipid metabolism, neural tissue function, and skin barrier formation.

[Learn More](#)

Parameter	Value
Product Name	FA2H Rabbit pAb
Remarks/Alias	FAAH, FAH1, SCS7, SPG35, FAXDC1, FA2H
Species	Human
Gene ID	79152 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 95-170 of human FA2H (NP_077282.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7L5A8
Protein Weight	43kDa
Shipping	Ice bag