

Anti-Faah1 Monoclonal Antibody For Wb, If/Icc, Elisa - 000519

Item Number: CM0012154

Introduction

Rabbit monoclonal FAAH1 antibody validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat FAAH1 (63 kDa). Ideal for endocannabinoid and lipid signaling research.

[Learn More](#)

Parameter	Value
Product Name	FAAH1 Rabbit mAb
Remarks/Alias	PSAB; FAAH-1; FAAH1
Species	Human
Gene ID (Human)	2166
Gene ID	2166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-579 of human FAAH1 (O00519).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00519
Protein Weight	63 kDa
Shipping	Ice bag