

Anti-Faah1 Rabbit Polyclonal Antibody For Wb, Elisa - 000519

Item Number: CM0015101

Introduction

High-quality rabbit polyclonal antibody against fatty-acid amide hydrolase 1 (FAAH1), validated for Western blot and ELISA. Reacts with human and mouse FAAH1. Ideal for endocannabinoid system research.

[Learn More](#)

Parameter	Value
Product Name	FAAH1 Rabbit pAb
Remarks/Alias	PSAB; FAAH-1; FAAH1
Species	Human
GeneID	2166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-579 of human FAAH11 (NP_001432.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00519
Protein Weight	63kDa
Shipping	Ice bag