

Anti-Faah2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Gmr7

Item Number: CM0024405

Introduction

FAAH2 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA applications. Detects human, mouse, rat FAAH2. Suitable for studies on lipid metabolism, endocannabinoid signaling, and amidase activity.

[Learn More](#)

Parameter	Value
Product Name	FAAH2 Rabbit pAb
Remarks/Alias	AMDD; FAAH2
Species	Human
Gene ID	158584
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 323-532 of human FAAH2 (NP_777572.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6GMR7
Protein Weight	58kDa
Shipping	Ice bag