

Anti-Fabp1 Rabbit Monoclonal Antibody For Wb - P07148

Item Number: CM0003570



Introduction

Rabbit monoclonal antibody targeting human FABP1 (L-FABP), suitable for WB, IF/ICC, IF-P, and ELISA. Detects FABP1 in human, mouse, and rat samples. Involved in intracellular lipid transport and cholesterol uptake.

[Learn More](#)

Parameter	Value
Product Name	FABP1 Rabbit mAb
Remarks/Alias	FABPL; L-FABP; FABP1
Species	Human
Gene ID (Human)	2168
Gene ID	2168
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FABP1 (P07148))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07148
Protein Weight	14kDa
Shipping	Ice bag