

Anti-Fabp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P07148

Item Number: CM0028745



Introduction

Rabbit polyclonal antibody against human FABP1 (L-FABP) for Western blot, IF/ICC, and ELISA detection. Cross-reacts with mouse and rat. Suitable for lipid metabolism and liver function studies.

[Learn More](#)

Parameter	Value
Product Name	FABP1 Rabbit pAb
Remarks/Alias	FABPL; L-FABP; FABP1
Species	Human
GeneID (Human)	2168
GeneID	2168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human FABP1 (NP_001434.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07148
Protein Weight	14kDa
Shipping	Ice bag