

Anti-Fabp12 Polyclonal Antibody For Wb, Elisa - A6Nfh5

Item Number: CM0025921

Introduction

Anti-FABP12 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human and rat FABP12, a lipid transport protein expressed in retinoblastoma. Ideal for lipid metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	FABP12 Rabbit pAb
Remarks/Alias	FABP12
Species	Human
GeneID (Human)	646486
GeneID	646486
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human FABP12 (NP_001098751.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	A6NFH5
Protein Weight	16kDa
Shipping	Ice bag