


Anti-Fabp4 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15090

Item Number: CM0014080



Introduction

Rabbit polyclonal antibody against human FABP4 (adipocyte fatty acid-binding protein, aP2). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for lipid transport and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	FABP4 Rabbit pAb
Remarks/Alias	aP2; ALBP; AFABP; A-FABP; HEL-S-104; FABP4
Species	Human
Gene ID (Human)	2167
Gene ID	2167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human FABP4 (NP_001433.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15090
Protein Weight	15 kDa
Shipping	Ice bag