


Anti-Fabp5 Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q01469

Item Number: CM0004725



Introduction

FABP5 Rabbit Monoclonal Antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA in Human, Mouse, Rat samples. Carrier for long-chain fatty acids, endocannabinoids; involved in inflammation and keratinocyte differentiation.

[Learn More](#)

Parameter	Value
Product Name	FABP5 Rabbit mAb
Remarks/Alias	EFABP; KFABP; E-FABP; PAFABP; PA-FABP
Species	Human
GeneID	2171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human FABP5 (NP_001435.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01469
Protein Weight	15kDa
Shipping	Ice bag