


Anti-Fabp5 Rabbit Polyclonal Antibody For Wb And Elisa - Q01469

Item Number: CM0015443



Introduction

Rabbit polyclonal antibody against human FABP5 (E-FABP), validated for WB and ELISA. Detects human, mouse, and rat FABP5; intracellular fatty acid transporter involved in inflammation and keratinocyte differentiation.

[Learn More](#)

Parameter	Specification
Product Name	FABP5 Rabbit pAb
Remarks/Alias	EFABP; KFABP; E-FABP; PAFABP; PA-FABP; FABP5
Species	Human
Gene ID (Human)	2171
Gene ID	2171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FABP5 (NP_001435.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01469
Protein Weight	15kDa
Shipping	Ice bag