


Anti-Fabp5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q01469

Item Number: CM0024740



Introduction

High-quality rabbit polyclonal antibody against human FABP5 (E-FABP), validated for WB, IHC-P, IF/ICC, and ELISA with human, mouse, and rat samples. Ideal for fatty acid transport, inflammation, and keratinocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	FABP5 Rabbit pAb
Remarks/Alias	EFABP; KFABP; E-FABP; PAFABP; PA-FABP; FABP5
Species	Human
GeneID (Human)	2171
GeneID	2171
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human FABP5 (NP_001435.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01469
Protein Weight	15kDa
Shipping	Ice bag