


Anti-Fabpi (I-Fabp) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P12104

Item Number: CM0007881



Introduction

High-quality rabbit monoclonal antibody against human FABPI (I-FABP), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for intestinal fatty acid transport and metabolic research.

[Learn More](#)

Parameter	Description
Product Name	FABPI Rabbit mAb
Remarks/Alias	FABPI; I-FABP
Species	Human
GeneID (Human)	2169
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human FABPI (P12104)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12104
Protein Weight	15kDa
Shipping	Ice bag