


Anti-Fabpi Polyclonal Antibody For Wb, Elisa - P12104

Item Number: CM0024084



Introduction

High-quality rabbit polyclonal antibody against human FABPI (I-FABP), validated for Western blot and ELISA applications. Detects human, mouse, and rat intestinal fatty acid-binding protein. Ideal for lipid transport research.

[Learn More](#)

Parameter	Value
Product Name	FABPI Rabbit pAb
Remarks/Alias	FABPI; I-FABP
Species	Human
Gene ID (Human)	2169
Gene ID	2169
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human FABPI (NP_000125.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12104
Protein Weight	15kDa
Shipping	Ice bag