


Anti-Factor IX / F9 Polyclonal Antibody For Wb, If-P, Elisa - P00740

Item Number: CM0022917



Introduction

Rabbit polyclonal antibody against human Factor IX (F9, FIX, Christmas factor). Validated for WB, IF-P, ELISA; cross-reacts with mouse, rat. Suitable for coagulation research and hemostasis diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Factor IX / F9 Rabbit pAb
Remarks/Alias	FIX; P19; PTC; HEMB; THPH8; F9 p22; Factor IX / F9
Species	Human
GeneID (Human)	2158
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 29-192 of human Factor IX / F9 (NP_000124.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00740
Protein Weight	52kDa
Shipping	Ice bag