


Anti-Fadd Polyclonal Antibody For Wb, Ip, If-P, Elisa - Q61160

Item Number: CM0021085



Introduction

Rabbit polyclonal antibody targeting FADD (Mort1/FADD), raised against mouse FADD (aa 1-205). Suitable for WB, IP, IF-P, ELISA; cross-reacts with human, mouse, rat. Key DISC adapter in apoptosis and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	FADD Rabbit pAb
Remarks/Alias	Mort1/FADD; FADD
Species	Mouse
GeneID (Human)	8772
GeneID	14082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of mouse FADD (NP_034305.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q61160
Protein Weight	23kDa
Shipping	Ice bag