


Anti-Fadd Rabbit Monoclonal Antibody For Wb, Elisa - Q13158

Item Number: CM0003117



Introduction

Rabbit monoclonal antibody targeting human FADD (GIG3/MORT1), validated for WB and ELISA applications. Ideal for apoptosis and innate immunity research. - Q13158

[Learn More](#)

Parameter	Value
Product Name	FADD Rabbit mAb
Remarks/Alias	GIG3; IMD90; MORT1; FADD
Species	Human
GeneID (Human)	8772
GeneID	8772
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-150 of human FADD (Q13158).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13158
Protein Weight	23kDa
Shipping	Ice bag