

## Anti-Fads1 Monoclonal Antibody For Wb And Elisa - O60427

Item Number: CM0002739

### Introduction

Anti-FADS1 rabbit monoclonal antibody for WB and ELISA, validated in human, mouse, and rat. Targets acyl-CoA desaturase involved in PUFA metabolism and inflammatory response.

[Learn More](#)

| Parameter        | Value                                                               |
|------------------|---------------------------------------------------------------------|
| Product Name     | FADS1 Rabbit mAb                                                    |
| Remarks/Alias    | D5D; TU12; FADS6; FADSD5; LLCDL1; FADS1                             |
| Species          | Human                                                               |
| GeneID           | 3992                                                                |
| Immunogen        | Synthetic peptide within amino acids 50-150 of human FADS1 (O60427) |
| Source           | Rabbit                                                              |
| Category         | Monoclonal Antibodies                                               |
| Application      | WB, ELISA                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                   |
| SWISS            | O60427                                                              |
| Protein Weight   | 52kDa                                                               |
| Shipping         | Ice bag                                                             |