

Anti-Fads3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y5Q0

Item Number: CM0027916

Introduction

FADS3 Rabbit pAb is a polyclonal antibody generated against recombinant FADS3 protein. It detects human, mouse, and rat FADS3, with applications in WB and ELISA. Ideal for studying fatty acid desaturase 3 in sphingolipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	FADS3 Rabbit pAb
Remarks/Alias	CYB5RP; LLCDL3; FADS3
Species	Human
GeneID (Human)	3995
GeneID	3995
Immunogen	Recombinant protein corresponding to amino acids 1-130 of human FADS3 (NP_068373.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5Q0
Protein Weight	51kDa
Shipping	Ice bag