


Anti-Faf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Unn5

Item Number: CM0028595



Introduction

CamelBio supplies an Anti-FAF1 rabbit polyclonal antibody for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat; detects a 74 kDa band. Useful for apoptosis, DNA replication, and nuclear protein research.

[Learn More](#)

Parameter	Value
Product Name	FAF1 Rabbit pAb
Remarks/Alias	hFAF1, CGI-03, HFAF1s, UBXD12, UBXN3A, FAF1
Species	Human
Gene ID (Human)	11124
Gene ID	11124
Immunogen	A synthetic peptide of human FAF1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNN5
Protein Weight	74kDa
Shipping	Ice bag