


## Anti-Faf2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Cs3

Item Number: CM0019608



### Introduction

FAF2 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, rat FAF2. Ideal for ERAD, lipid droplet, and stress granule research.

[Learn More](#)

| Parameter        | Value                     |
|------------------|---------------------------|
| Product Name     | FAF2 Rabbit pAb           |
| Remarks/Alias    | ETEa; UBxD8; UBxN3B; FAF2 |
| Species          | Human                     |
| Gene ID (Human)  | 23197                     |
| Gene ID          | 23197                     |
| Immunogen        | Recombinant protein       |
| Source           | Rabbit                    |
| Category         | Polyclonal Antibodies     |
| Application      | WB, IF/ICC, ELISA         |
| Cross Reactivity | Human, Mouse, Rat         |
| SWISS            | Q96CS3                    |
| Protein Weight   | 53kDa                     |
| Shipping         | Ice bag                   |