

Anti-Fah Polyclonal Antibody For Wb, If/Icc, Elisa - P16930

Item Number: CM0025477

Introduction

Anti-FAH polyclonal antibody (Rabbit) validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets fumarylacetoacetase, a key enzyme in tyrosine catabolism. Ideal for metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	FAH Rabbit pAb
Remarks/Alias	FAH
Species	Human
Gene ID	2184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human FAH (NP_000128.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16930
Protein Weight	46kDa
Shipping	Ice bag